



Agilent's New, advanced  
**ZORBAX Rapid Resolution HT (RRHT)**  
HPLC Columns

*The ideal combination of*  
**Speed + Resolution + Sensitivity + Stability**



**600 Bar**



---

---

***Thoroughbred performance  
saves up to 95% in HPLC analysis time***

---

---

**LCGC**

*Exclusive solutions for Innovative minds*



**Agilent Technologies**

— Authorised Distributor

# Agilent's Advanced High Pressure ZORBAX RRHT 1.8 $\mu$ m LC Columns:

## Helping you overcome your most daunting separation challenges.

You asked for choice in column dimension, ID, and bonded phase.

We offer an unsurpassed variety of columns to accommodate virtually any application.

Choose from the following configurations optimized for performance, reliability, ruggedness, and long life:

### ZORBAX RRHT 1.8 $\mu$ m LC columns for high-pressure use (maximum pressure: 600 bar, 9000 PSI)

Best for fast and high-resolution LC.

	SB-C18	SB-C8	Eclipse XDB-C-18	Eclipse XDB-C8	Extend-C18	SB-CN
4.6 x 150 mm	829975-902					
4.6 x 100 mm	828975-902	828975-906	928975-902		728975-902	
4.6 x 50 mm	827975-902	827975-906	927975-902	927975-906	727975-902	827975-905
4.6 x 30 mm	824975-902	824975-906	924975-902	924975-906	724975-902	824975-905
4.6 x 20 mm	826975-902	826975-906	926975-902	926975-906	726975-902	
3.0 x 150 mm	829975-302					
3.0 x 100 mm	828975-302	828975-306	928975-302		728975-302	
3.0 x 50 mm	827975-302	827975-306	927975-302	927975-306	727975-302	827975-305
3.0 x 30 mm	824975-302	824975-306	924975-302	924975-306	724975-302	824975-305
3.0 x 20 mm	826975-302	826975-306	926975-302	926975-306	726975-302	
2.1 x 150 mm	820700-902	820700-905				820700-905
2.1 x 100 mm	828700-902	828700-906	928700-902	928700-906	728700-902	828700-905
2.1 x 50 mm	827700-902	827700-906	927700-902	927700-906	727700-902	827700-905
2.1 x 30 mm	824700-902	824700-906	924700-902	924700-906	724700-902	824700-905
2.1 x 20 mm	826700-902	826700-906	926700-902	926700-906	726700-902	
1.0 x 100 mm	828600-902	828600-906	928600-906	928600-902	728600-902	828600-905
1.0 x 50 mm	822600-902	822600-906	922600-902	922600-906	722600-902	822600-905

### ZORBAX RRHT 1.8 $\mu$ m columns and cartridges (maximum pressure: 400 bar, 6000 PSI)

Best for fast LC.

	Eclipse XDB-18	Eclipse XDB-C18 3pk	Eclipse XDB-C8	SB-C18	SB-C18 3pk	SB-C8	Extend-C18
Standard Columns							
4.6 x 50 mm	922975-902	922975-932	922975-906	822975-902	822975-932	822975-906	722975-902
2.1 x 50 mm	922700-902	922700-932		822700-902	822700-932		
Cartridge Columns (require hardware kit 820555-901)							
4.6 x 50 mm	925975-902	925975-932		825975-902	825975-932		
4.6 x 30 mm	923975-902	923975-932		823975-902	823975-932		
4.6 x 15 mm	921975-902	921975-932		821975-902	821975-932		
2.1 x 50 mm	925700-902	925700-932		825700-902	825700-932		
2.1 x 30 mm	923700-902	923700-932		823700-902	823700-932		
2.1 x 15 mm	921700-902	921700-932		821700-902	821700-932		



## LCGC CHROMATOGRAPHY SOLUTIONS PVT LTD

E-mail : [info@lccgindia.com](mailto:info@lccgindia.com) Web : [www.lccgindia.com](http://www.lccgindia.com)

**Head Office** **Hyderabad**  
 LCGC Towers, 161/166, Paigah Colony,  
 Near Anand Theatre, S. P. Road, Secunderabad - 500 003.  
 Phone : 040-66315242, 66315243  
 Reliance Mobile : 9347015242  
 Fax : 040-55329000  
 E-mail : [INFO@LCGCINDIA.COM](mailto:INFO@LCGCINDIA.COM)

**Main Branch** **Mumbai**  
 104, Champaklal Industrial Estate,  
 Near Sion Telephone Exchange, Sion (East), Mumbai- 400 022.  
 Phone : 022-24073068, 24091124  
 Reliance Mobile : 9323415242  
 Fax : 022-24074562  
 E-mail : [MUMBAI@LCGCINDIA.COM](mailto:MUMBAI@LCGCINDIA.COM)

**Branches** **Ahmedabad** M. P. Sriram : 98253 17798  
 B-3/10,3rd Floor, Shanti Apartment,  
 Near Ashish School, Bodakdev,  
 Ahmedabad-380 051.  
 Phone : 079- 26871043  
 Reliance Mobile : 9374015242  
 Fax : 079-2687 2926  
 E-mail : [AHMEDABAD@LCGCINDIA.COM](mailto:AHMEDABAD@LCGCINDIA.COM)

**Branches** **Bangalore** M.N. Suresh : 98440 21071  
 MS-1, A-Block, Flat No-1, Nandini Layout,  
 Bangalore-560 096.  
 Phone : 080-23595827  
 Reliance Mobile : 9342855242  
 Fax : 080-23595827  
 E-mail : [BANGALORE@LCGCINDIA.COM](mailto:BANGALORE@LCGCINDIA.COM)

**Chennai** Rajesh Krishnan : 98410 56364  
 18, 2nd Main Road, Krishna Nagar, Chrompet,  
 Chennai - 600 044.  
 Phone : 044-22238259  
 Fax : 044-22233483  
 E-mail : [CHENNAI@LCGCINDIA.COM](mailto:CHENNAI@LCGCINDIA.COM)

**New Delhi** C.K.N.Rajan : 98103 87404  
 S-5, Kiran Jyothi Complex, Nr. Amar Jyothi  
 Housing Society, Mayur Vihar, Phase-1,  
 New Delhi-110 091.  
 Phone : 011-22756195  
 Reliance Mobile : 9313415242  
 Fax : 011-22794814  
 E-mail : [DELHI@LCGCINDIA.COM](mailto:DELHI@LCGCINDIA.COM)

**Branches** **Chandigarh (Baddi)** Pankaj Rai.  
 H.No.3199/2, MIG Cat-II, Sector-44D, Chandigarh-160 047 (U.T.)  
 Reliance Mobile : 9357448808  
 E-mail : [PANKAJ@LCGCINDIA.COM](mailto:PANKAJ@LCGCINDIA.COM)

**Pune** Purshottam Pawar - 9373015242  
 E-18, Rahul Complex, Nr. Krishna Hospital,  
 Paud Road, Kothrud, Pune-411 038.  
 Telefax : 020-25468838  
 E-mail : [PUNE@LCGCINDIA.COM](mailto:PUNE@LCGCINDIA.COM)

**Ankleshwar (Surat)** Ketan R. Patel: 9824149149  
 12-B, Aataman Park Society, Near Haridarshan Row House,  
 New Adajan-Pal Road, Surat-395 003. Gujarat  
 E-mail : [KETAN@LCGCINDIA.COM](mailto:KETAN@LCGCINDIA.COM)

**Kolkatta** **Tanmoy** Mitra : 94330 85262  
 Phone : 033-22271107 / 1108  
 E-mail : [TANMOY@CHROMSPEC.NET](mailto:TANMOY@CHROMSPEC.NET)  
 Fax : 033-22271108

**Indore** Ashish Dayma : 98260 84686  
 101 Gurukripa Apartment, 25 Kalpana lok Colony,  
 Behind Anand Bazar Area, Indore-452 018. (M.P.)  
 Telefax : 0731-2561037  
 E-mail : [ASHISH@LCGCINDIA.COM](mailto:ASHISH@LCGCINDIA.COM)